

**Peristaltic pump tubing, Tygon LMT-55 with 2 colour-coded bridges****Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges**

With 2 colour-coded bridges (international standards). Shore: 55. Fits most common peristaltic pump models, e.g. Verder, Ismatec, Watson-Marlow, IKA and Heidolph. Standard tubing, transparent. Length 380 mm. Long working lifespan. Good resistance to acids and alkalis. Food-safe.



Code	Description	Packaging
LLG09208202	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-red, Int. diam. 0.19 mm, Wall thickness 0.9 mm	12 pz.
LLG09208207	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour black-black, Int. diam. 0.76 mm, Wall thickness 0.85 mm	12 pz
LLG09208222	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour black-white, Int. diam. 3.18 mm, Wall thickness 0.85 mm	12 pz
LLG09208213	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour blue-blue, Int. diam. 1.65 mm, Wall thickness 0.85 mm	12 pz
LLG09208212	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour yellow-yellow, Int. diam. 1.42 mm, Wall thickness 0.85 mm	12 pz
LLG09208214	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour green-green, Int. diam. 1.85 mm, Wall thickness 0.85 mm	12 pz
LLG09208211	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour grey-grey, Int. diam. 1.3 mm, Wall thickness 0.85 mm	12 pz
LLG09208215	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour purple-purple, Int. diam. 2.06 mm, Wall thickness 0.85 mm	12 pz
LLG09208201	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-black, Int. diam. 0.13 mm, Wall thickness 0.9 mm	12 pz.
LLG09208203	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-blue, Int. diam. 0.25 mm, Wall thickness 0.9 mm	12 pz.
LLG09208204	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-green, Int. diam. 0.38 mm, Wall thickness 0.9 mm	12 pz.
LLG09208208	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-orange, Int. diam. 0.89 mm, Wall thickness 0.85 mm	12 pz
LLG09208206	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-white, Int. diam. 0.64 mm, Wall thickness 0.9 mm	12 pz
LLG09208205	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour orange-yellow, Int. diam. 0.51 mm, Wall thickness 0.9 mm	12 pz



Code	Description	Packaging
LLG09208216	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour purple-black, Int. diam. 2.29 mm, Wall thickness 0.85 mm	12 pz
LLG09208217	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour purple-orange, Int. diam. 2.54 mm, Wall thickness 0.85 mm	12 pz
LLG09208218	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour purple-white, Int. diam. 2.79 mm, Wall thickness 0.85 mm	12 pz
LLG09208210	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour red-red, Int. diam. 1.14 mm, Wall thickness 0.85 mm	12 pz
LLG09208209	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour white-white, Int. diam. 1.02 mm, Wall thickness 0.85 mm	12 pz
LLG09208219	Peristaltic pump tubing, Tygon® LMT-55 with 2 colour-coded bridges, Colour yellow-blue, Int. diam. 1.52 mm, Wall thickness 0.85 mm	12 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

