

**Vampire Sampler****Bürkle****Vampire Sampler**

The self-priming peristaltic pump is also suitable for dirty or particle-containing substances with a viscosity of max. 500 mPas. Ideal for sampling from barrels, tanks, wells, water separators, waste water and sewage treatment plants. The samples are fed directly into the sample container via the closed hose system.

- Battery-powered hose pump
- Easy cleaning and exchangeable hoses prevent impurities
- Suction depth up to 5 metres

Scope of supply: Drive unit, pump head, pump hose (30 cm, Marprene), suction hose (5 m, PE), hose connector (PP), sample bottle (250 ml, LDPE), bottle basket, two rechargeable batteries, battery charger



Code	Description	Packaging
LLG06284634	Accessories for Vampire Sampler, Type Bottle holder for 500 ml bottle	1 pz.
LLG06233558	Accessories for Vampire Sampler, Type Extension tubing, PE, 10 m	10 m
LLG09303851	Vampire Sampler, Type Vampire sampler, complete	1 pz.
LLG09303854	Wide-necked bottle, LDPE, transparent, Capacity 250 ml	1 pz.
LLG09303855	Wide-necked bottle, LDPE, transparent, Capacity 500 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

