

**Tubes, glass, DURAN, culture, screw cap****DWK Life Sciences****Tubes, glass, DURAN®, culture, screw cap**

DIN thread. Complete with red PBT screw cap with PTFE-faced sealing wad. Thus, it is suitable for the cultivation of microorganisms. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable.



Code	Description	Packaging
LLG09010012	Tubes, glass, DURAN®, culture, screw cap, Diam. 12 mm, Height 100 mm, Neck thread 14 GL	50 pz
LLG09010013	Tubes, glass, DURAN®, culture, screw cap, Diam. 13 mm, Height 100 mm, Neck thread 14 GL	50 pz
LLG09010014	Tubes, glass, DURAN®, culture, screw cap, Diam. 16 mm, Height 150 mm, Neck thread 18 GL	50 pz
LLG09010016	Tubes, glass, DURAN®, culture, screw cap, Diam. 16 mm, Height 160 mm, Neck thread 18 GL	50 pz
LLG09010018	Tubes, glass, DURAN®, culture, screw cap, Diam. 18 mm, Height 180 mm, Neck thread 18 GL	50 pz
LLG09010019	Tubes, glass, DURAN®, culture, screw cap, Diam. 20 mm, Height 150 mm, Neck thread 18 GL	50 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

