

**Flex tubing, PFA****BOLA****Flex tubing, PFA**

Corrugated tubing with circular corrugations around the longitudinal axis and cylindrical tubing ends with a length of 40mm which can for example be connected directly to fittings, stopcocks or hose connectors. Together with BOLA Laboratory Screw Joints, the connection is tight and suitable for vacuum.

- Highly flexible
- Tight bending radius only causes little cross-section reduction
- Non-porous
- Transparent
- Ideal for connections under vibrations
- Usable with a small bending radius
- For compensation of thermal expansions
- For easy handling of liquids



Code	Description	Packaging
LLG09205713	Flex tubing, PFA, Int. diam. 10 mm, Ext. diam. 12 mm, Wall thickness 1 mm, Length 0.5 m, Size 13 NW	1 pz.
LLG09205733	Flex tubing, PFA, Int. diam. 10 mm, Ext. diam. 12 mm, Wall thickness 1 mm, Length 1 m, Size 13 NW	1 pz.
LLG09205714	Flex tubing, PFA, Int. diam. 12 mm, Ext. diam. 14 mm, Wall thickness 1 mm, Length 0.5 m, Size 14 NW	1 pz.
LLG09205719	Flex tubing, PFA, Int. diam. 20.9 mm, Ext. diam. 25.4 mm, Wall thickness 2.25 mm, Length 0.5 m, Size 23 NW	
LLG09205734	Flex tubing, PFA, Int. diam. 12 mm, Ext. diam. 14 mm, Wall thickness 1 mm, Length 1 m, Size 14 NW	1 pz.
LLG09205739	Flex tubing, PFA, Int. diam. 20.9 mm, Ext. diam. 25.4 mm, Wall thickness 2.25 mm, Length 1 m, Size 23 NW	
LLG09205712	Flex tubing, PFA, Int. diam. 8 mm, Ext. diam. 10 mm, Wall thickness 1 mm, Length 0.5 m, Size 10 NW	1 pz.
LLG09205732	Flex tubing, PFA, Int. diam. 8 mm, Ext. diam. 10 mm, Wall thickness 1 mm, Length 1 m, Size 10 NW	1 pz.
LLG09205741	Flex tubing, PFA, Int. diam. 8 mm, Ext. diam. 10 mm, Wall thickness 1 mm, Length 0.25 m, Size 10 NW	1 pz.
LLG09205717	Flex tubing, PFA, Int. diam. 16 mm, Ext. diam. 18 mm, Wall thickness 1 mm, Length 0.5 m, Size 19 NW	
LLG09205737	Flex tubing, PFA, Int. diam. 16 mm, Ext. diam. 18 mm, Wall thickness 1 mm, Length 1 m, Size 19 NW	
LLG09205715	Flex tubing, PFA, Int. diam. 14 mm, Ext. diam. 16 mm, Wall thickness 1 mm, Length 0.5 m, Size 16 NW	
LLG09205735	Flex tubing, PFA, Int. diam. 14 mm, Ext. diam. 16 mm, Wall thickness 1 mm, Length 1 m, Size 16 NW	1 pz.
LLG09205718	Flex tubing, PFA, Int. diam. 17.5 mm, Ext. diam. 20 mm, Wall thickness 1.25 mm, Length 0.5 m, Size 21 NW	
LLG09205738	Flex tubing, PFA, Int. diam. 17.5 mm, Ext. diam. 20 mm, Wall thickness 1.25 mm, Length 1 m, Size 21 NW	
LLG09205731	Flex tubing, PFA, Int. diam. 6 mm, Ext. diam. 8 mm, Wall thickness 1 mm, Length 1 m, Size 8 NW	1 pz.
LLG06237188	Flex tubing, PFA, Int. diam. 8 mm, Ext. diam. 10 mm, Wall thickness 1 mm, Length 2.5 m, Size 10 NW	1 pz.
LLG09205740	Flex tubing, PFA, Int. diam. 10 mm, Ext. diam. 12 mm, Wall thickness 1 mm, Length 2.5 m, Size 13 NW	1 pz.
LLG09205710	Flex tubing, PFA, Int. diam. 2 mm, Ext. diam. 4 mm, Wall thickness 1 mm, Length 0.5 m, Size 4.5 NW	



Code	Description	Packaging
LLG09205730	Flex tubing, PFA, Int. diam. 2 mm, Ext. diam. 4 mm, Wall thickness 1 mm, Length 1 m, Size 4.5 NW	
LLG09205711	Flex tubing, PFA, Int. diam. 6 mm, Ext. diam. 8 mm, Wall thickness 1 mm, Length 0.5 m, Size 8 NW	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

