

**GL stopcocks, PTFE****BOLA****GL stopcocks, PTFE**

Stopcocks with cylindrical taps and GL screwthread made of pure PTFE, giving almost total chemical resistance. Can be used at working pressures of up to 6 bar when used in combination with Bola laboratory screwthread connections (at room temperature).

A further advantage is that they can be used as expansion/reduction adapters to connect tubing of different diameters. For easy tubing connection we recommend accessory Bola laboratory screwthread connectors - please order separately.



Code	Description	Packaging
LLG09207153	GL stopcocks, PTFE, 2-way, Bore diam. 4 mm, Neck thread 14 GL	1 pz.
LLG09207156	Stopcocks, 3-way, 4mm boreGL 14, PTFE	1 pz.
LLG09207154	GL stopcocks, PTFE, 2-way, Bore diam. 6 mm, Neck thread 18 GL	1 pz.
LLG09207157	Stopcocks, 3-way, 6mm boreGL 18, PTFE	1 pz.
LLG09207158	Stopcocks, 3-way, 6mm boreGL 25, PTFE	1 pz.
LLG09207155	GL stopcocks, PTFE, 2-way, Bore diam. 8 mm, Neck thread 25 GL	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

