

**Microlitre syringes, 800 series, with cemented (N) or removable needles (RN)**

Point Style 2 (pst2): 12° bevelled non-coring needle point recommended for septum penetration. Ideal for gas chromatographic applications.

Point Style 3 (pst3): Blunt needle point (90°) for use with HPLC injection valves. Also recommended for applications where exact dosing is required (e.g. thin layer chromatographs).

Point Style 4 (pst4): 10-12° bevelled needle point recommended for life science applications; sharp point for animal injection.

Point Style 5 (pst5): Conical needle with side port for penetration of septa, thin-gauged vinyls and plastics without coring; minimizes septum damage

Point Style AS (pstAS): Special conical style needle point (8° taper) designed to withstand the demands of multiple injections; exclusively used on GC autosampler syringes

With metal handle. Plungers are individually fitted, therefore cannot be interchanged and are not available as replacement parts. Cemented needles (N) or removable needle (RN).



Code	Description	Packaging
LLG06222013	Microlitre syringes, Type 801 RN*, 10 µl, Gauge, Point style, Needle length mm	1 pz
LLG09221110	Microlitre syringes, Type 801 N, 10 µl, Gauge 26s, Point style 2, Needle length 51 mm	1 pz
LLG09221116	Microlitre syringes, Type 801 RN, 10 µl, Gauge 26s, Point style 2, Needle length 51 mm	1 pz
LLG06231153	Microlitre syringes, Type 810 RN*, 100 µl, Gauge, Point style, Needle length mm	1 pz
LLG09221113	Microlitre syringes, Type 810 N, 100 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221119	Microlitre syringes, Type 810 RN, 100 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221111	Microlitre syringes, Type 802 N, 25 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG06232019	Microlitre syringes, Type 85 RN B/P, 5 µl, Gauge, Point style, Needle length mm	1 pz.
LLG06803384	Microlitre syringes, Type 85 RN*, 5 µl, Gauge, Point style, Needle length mm	1 pz
LLG09221115	Microlitre syringes, Type 85 RN, 5 µl, Gauge 26s, Point style 2, Needle length 51 mm	1 pz
LLG09221112	Microlitre syringes, Type 805 N, 50 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221118	Microlitre syringes, Type 805 RN, 50 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221172	Microlitre syringes, Type 801 RN B/P, 10 µl, Gauge, Point style, Needle length mm	1 pz
LLG09221117	Microlitre syringes, Type 802 RN, 25 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221120	Microlitre syringes, Type 825 RN, 250 µl, Gauge 22s, Point style 2, Needle length 51 mm	1 pz
LLG09221105	Microlitre syringes, Type 85 N, 5 µl, Gauge 26s, Point style 2, Needle length 51 mm	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

