

**Needles for LT / TLL / TLLX syringes****Hamilton****Needles for LT / TLL / TLLX syringes**

Point Style 2 (pst2): 12° bevelled non-coring needle point recommended for septum penetration. Ideal for gas chromatographic applications.

Point Style 3 (pst3): Blunt needle point (90°) for use with HPLC injection valves. Also recommended for applications where exact dosing is required (e.g. thin layer chromatographs).

Point Style 4 (pst4): 10-12° bevelled needle point recommended for life science applications; sharp point for animal injection.

Point Style 5 (pst5): Conical needle with side port for penetration of septa, thin-gauged vinyls and plastics without coring; minimizes septum damage

Point Style AS (pstAS): Special conical style needle point (8° taper) designed to withstand the demands of multiple injections; exclusively used on GC autosampler syringes

Removable. With Luer tip made of Kel-F. All hypodermic needles available according to specification in 10 mm to max. 1000 mm lengths. Further point styles and diameters available on request. Needle length: 51 mm



Code	Description	Packaging
LLG09221726	Needles for LT / TLL / TLLX syringes, Type KF 726 pst3, Ext. diam. 0.46 mm, Int. diam. 0.26 mm, Gauge 26, Point style 3 (HPLC)	6 pz
LLG09221609	Needles for LT / TLL / TLLX syringes, Type KF 726s pst2, Ext. diam. 0.47 mm, Int. diam. 0.13 mm, Gauge 26s, Point style 2 (GC)	6 pz
LLG09221709	Needles for LT / TLL / TLLX syringes, Type KF 726s pst3, Ext. diam. 0.47 mm, Int. diam. 0.13 mm, Gauge 26s, Point style 3 (HPLC)	6 pz
LLG09221625	Needles for LT / TLL / TLLX syringes, Type KF 725 pst2, Ext. diam. 0.51 mm, Int. diam. 0.26 mm, Gauge 25s, Point style 2 (GC)	6 pz
LLG09221735	Needles for LT / TLL / TLLX syringes, Type KF 725 pst3, Ext. diam. 0.51 mm, Int. diam. 0.26 mm, Gauge 25, Point style 3 (HPLC)	6 pz
LLG09221624	Needles for LT / TLL / TLLX syringes, Type KF 724 pst2, Ext. diam. 0.57 mm, Int. diam. 0.31 mm, Gauge 24s, Point style 2 (GC)	6 pz
LLG09221724	Needles for LT / TLL / TLLX syringes, Type KF 724 pst3, Ext. diam. 0.57 mm, Int. diam. 0.31 mm, Gauge 24, Point style 3 (HPLC)	6 pz
LLG09221607	Needles for LT / TLL / TLLX syringes, Type KF 722s pst2, Ext. diam. 0.72 mm, Int. diam. 0.15 mm, Gauge 22s, Point style 2 (GC)	6 pz
LLG09221622	Needles for LT / TLL / TLLX syringes, Type KF 722 pst2, Ext. diam. 0.72 mm, Int. diam. 0.41 mm, Gauge 22, Point style 2 (GC)	6 pz
LLG09221707	Needles for LT / TLL / TLLX syringes, Type KF 722s pst3, Ext. diam. 0.72 mm, Int. diam. 0.15 mm, Gauge 22s, Point style 3 (HPLC)	6 pz
LLG09221623	Needles for LT / TLL / TLLX syringes, Type KF 723 pst2, Ext. diam. 0.64 mm, Int. diam. 0.34 mm, Gauge 23s, Point style 2 (GC)	6 pz
LLG09221723	Needles for LT / TLL / TLLX syringes, Type KF 723 pst3, Ext. diam. 0.64 mm, Int. diam. 0.34 mm, Gauge 23, Point style 3 (HPLC)	6 pz
LLG09221720	Needles for LT / TLL / TLLX syringes, Type KF 720 pst3, Ext. diam. 0.91 mm, Int. diam. 0.6 mm, Gauge 20, Point style 3 (HPLC)	6 pz
LLG09221621	Needles for LT / TLL / TLLX syringes, Type KF 721 pst2, Ext. diam. 0.82 mm, Int. diam. 0.51 mm, Gauge 21, Point style 2 (GC)	6 pz



Code	Description	Packaging
LLG09221721	Needles for LT / TLL / TLLX syringes, Type KF 721 pst3, Ext. diam. 0.82 mm, Int. diam. 0.51 mm, Gauge 21, Point style 3 (HPLC)	6 pz
LLG09221620	Needles for LT / TLL / TLLX syringes, Type KF 720 pst2, Ext. diam. 0.91 mm, Int. diam. 0.6 mm, Gauge 20, Point style 2 (GC)	6 pz
LLG09221722	Needles for LT / TLL / TLLX syringes, Type KF 722 pst3, Ext. diam. 0.72 mm, Int. diam. 0.41 mm, Gauge 22, Point style 3 (HPLC)	6 pz
LLG09221626	Needles for LT / TLL / TLLX syringes, Type KF 726 pst2, Ext. diam. 0.46 mm, Int. diam. 0.26 mm, Gauge 26, Point style 2 (GC)	6 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

