



Tubing, PTFE

Tubing, PTFE

Highly chemical resistant and withstands temperatures from -200°C to 250°C.



Code	Description	Packaging
LLG09205609	Tubing, PTFE, Int. diam. 0.5 mm, Ext. diam. 1.6 mm, Wall thickness 0.55 mm	1 pz.
LLG09205608	Tubing, PTFE, Int. diam. 0.5 mm, Ext. diam. 1 mm, Wall thickness 0.25 mm	1 pz.
LLG09205680	Tubing, PTFE, Int. diam. 8 mm, Ext. diam. 10 mm, Wall thickness 1 mm	1 pz.
LLG09205681	Tubing, PTFE, Int. diam. 8 mm, Ext. diam. 9.5 mm, Wall thickness 0.75 mm	
LLG09205682	Tubing, PTFE, Int. diam. 8 mm, Ext. diam. 11 mm, Wall thickness 1.5 mm	
LLG09205685	Tubing, PTFE, Int. diam. 9 mm, Ext. diam. 11 mm, Wall thickness 1 mm	
LLG09205684	Tubing, PTFE, Int. diam. 8 mm, Ext. diam. 12 mm, Wall thickness 1.6 mm	1 m
LLG09205691	Tubing, PTFE, Int. diam. 10 mm, Ext. diam. 12 mm, Wall thickness 1 mm	1 pz.
LLG09205692	Tubing, PTFE, Int. diam. 10 mm, Ext. diam. 14 mm, Wall thickness 2 mm	
LLG09205616	Tubing, PTFE, Int. diam. 0.8 mm, Ext. diam. 1.6 mm, Wall thickness 0.4 mm	1 pz.
LLG09205619	Tubing, PTFE, Int. diam. 1 mm, Ext. diam. 1.6 mm, Wall thickness 0.4 mm	1 pz.
LLG09205621	Tubing, PTFE, Int. diam. 1 mm, Ext. diam. 2 mm, Wall thickness 0.5 mm	1 pz.
LLG09205630	Tubing, PTFE, Int. diam. 1.5 mm, Ext. diam. 2.1 mm, Wall thickness 0.3 mm	1 pz.
LLG06070190	Tubing, PTFE, Int. diam. 1 mm, Ext. diam. 3 mm, Wall thickness 1 mm	1 pz.
LLG09205644	Tubing, PTFE, Int. diam. 2 mm, Ext. diam. 3 mm, Wall thickness 0.5 mm	1 pz.
LLG09205634	Tubing, PTFE, Int. diam. 1.6 mm, Ext. diam. 3.2 mm, Wall thickness 0.8 mm	1 pz.
LLG09205633	Tubing, PTFE, Int. diam. 1.5 mm, Ext. diam. 3.5 mm, Wall thickness 1 mm	1 pz.
LLG09205645	Tubing, PTFE, Int. diam. 2 mm, Ext. diam. 4 mm, Wall thickness 1 mm	1 pz.
LLG09205657	Tubing, PTFE, Int. diam. 3 mm, Ext. diam. 4 mm, Wall thickness 0.5 mm	1 pz.
LLG09205658	Tubing, PTFE, Int. diam. 3 mm, Ext. diam. 5 mm, Wall thickness 1 mm	1 pz.



Code	Description	Packaging
LLG09205663	Tubing, PTFE, Int. diam. 4 mm, Ext. diam. 5 mm, Wall thickness 0.5 mm	1 pz.
LLG09205665	Tubing, PTFE, Int. diam. 4 mm, Ext. diam. 6 mm, Wall thickness 1 mm	1 pz.
LLG09205668	Tubing, PTFE, Int. diam. 5 mm, Ext. diam. 6 mm, Wall thickness 0.5 mm	1 pz.
LLG09205669	Tubing, PTFE, Int. diam. 5 mm, Ext. diam. 7 mm, Wall thickness 1 mm	
LLG09205672	Tubing, PTFE, Int. diam. 6 mm, Ext. diam. 8 mm, Wall thickness 1 mm	1 pz.
LLG09205678	Tubing, PTFE, Int. diam. 7 mm, Ext. diam. 9 mm, Wall thickness 1 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

