

**Culture tubes, glass DURAN, small, rimless****DWK Life Sciences****Culture tubes, glass DURAN®, small, rimless**

To DIN 38411. Rimless. Suitable for use with Kapsenberg caps. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable. Please order cap separately.



Code	Description	Packaging
LLG09010046	Culture tubes, glass DURAN®, small, rimless, Diam. 16 mm, Height 160 mm, Wall thickness 1.0 to 1.2 mm	100 pz
LLG09010048	Culture tubes, glass DURAN®, small, rimless, Diam. 18 mm, Height 180 mm, Wall thickness 1.0 to 1.2 mm	100 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

