

**Culture media bottles DURAN, glass, cylindrical****DWK Life Sciences****Culture media bottles DURAN®, glass, cylindrical**

Rimless. For use with Kapsenberg caps. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable.



Code	Description	Packaging
LLG09010060	Culture media bottles DURAN®, glass, cylindrical, 100 ml, Diam. 40 mm, Neck diam. 18 mm, Height 150 mm	1 pz
LLG09010070	Culture media bottles DURAN®, glass, cylindrical, 200 ml, Diam. 50 mm, Neck diam. 18 mm, Height 175 mm	1 pz
LLG09010050	Culture media bottles DURAN®, glass, cylindrical, 50 ml, Diam. 40 mm, Neck diam. 18 mm, Height 107 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

