

**Draining Rack, ABS****Bel-Art Products****Draining Rack, ABS**

Laboratory space is always in high demand. Scienceware® Lab-Aire® II Drying Racks are a practical solution for drying and storing fragile, frequently used labware. Benchtop Dryers lets you place your dryer wherever it will be most convenient and a sturdy base provides excellent balance and support. Wallmount dryers keep labware off the bench and out of the way as it dries. Easily installs on walls or standard pegboard using moulded-in tabs with key hole openings. Hardware not included.

- Rugged ABS plastic will never corrode and easily wipes clean
- Drip/storage trays have scalloped wells to hold long thin labware
- 3 different sizes of pegs can be moved and repositioned to hold a variety of labware including beakers, flasks, bottles, test tubes and cylinders
- With drain



Code	Description	Packaging
LLG09003301	Draining Rack, ABS, Width 375 mm, Depth 127 mm, Height 240* mm, Rods 76 mm 6 Quantity, Rods 101 mm 7 Quantity, Rods 127 mm 6 Quantity	1 pz.
LLG09003302	Draining Rack, ABS, Width 375 mm, Depth 178 mm, Height 390 mm, Rods 76 mm 12 Quantity, Rods 101 mm 14 Quantity, Rods 127 mm 12 Quantity	1 pz.
LLG09003303	Draining Rack, ABS, Width 375 mm, Depth 127 mm, Height 594* mm, Rods 76 mm 18 Quantity, Rods 101 mm 21 Quantity, Rods 127 mm 18 Quantity	1 pz.
LLG09003304	Draining Rack, ABS, Width 375 mm, Depth 127 mm, Height 417* mm, Rods 76 mm 12 Quantity, Rods 101 mm 14 Quantity, Rods 127 mm 12 Quantity	1 pz.
LLG09003306	Draining Rack, ABS, Width 375 mm, Depth 178 mm, Height 570 mm, Rods 76 mm 18 Quantity, Rods 101 mm 21 Quantity, Rods 127 mm 18 Quantity	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

