

**96-Well PCR Rack, PP**

Heathrow Scientific

**96-Well PCR Rack, PP**

Removable hinged lid lies flat. Perfect for pre- or post-PCR sampling. Sturdy polypropylene rack accommodates individual 0.2ml tubes or strips of 8 or 12 tubes. Tube wells are easily identified with embossed numbers and letters. Dimensions (L x W x H): 130 x 98 x 33mm. Autoclavable.



Code	Description	Packaging
LLG06256114	96-Well PCR Rack, PP, For tubes 0.2 ml, Colour green	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

