

**Spray bottles LaboPlast, PE / PP****Bürkle****Spray bottles LaboPlast®, PE / PP**

Chemical resistant to a large extent. Trigger-action spray pump with continuously adjustable spray nozzle.

- Parts in contact with media made of PP, PE and stainless steel
- Nozzle diameter 0.6mm; nozzle is continuously adjustable
- Amount sprayed per stroke: 1.2ml \pm 0.1ml



Code	Description	Packaging
LLG09223160	Spray bottles LaboPlast®, PE / PP, 1000 ml, Ext. diam. 98 mm, Neck diam. int. 17.5 mm, Height 295 mm, Weight 111 g	1 pz.
LLG09223152	Spray bottles LaboPlast®, PE / PP, 250 ml, Ext. diam. 60 mm, Neck diam. int. 18 mm, Height 220 mm, Weight 54 g	1 pz.
LLG09223155	Spray bottles LaboPlast®, PE / PP, 500 ml, Ext. diam. 75 mm, Neck diam. int. 18 mm, Height 240 mm, Weight 79 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

