



Wash bottles Nalgene Type 2403, FEP, with screw cap, ETFE

Thermo Scientific

Wash bottles Nalgeneä Type 2403, FEP, with screw cap, ETFE

Screw cap, jet and feed tube are one-piece moulded from ETFE. Ideal for contamination-free work and suitable for trace element analysis.



Code	Description	Packaging
LLG09223234	Wash bottles Nalgene Type 2403, FEP, with screw cap, ETFE, 1000 ml, Screw cap diam. 38-430 mm	1 pz.
LLG09223231	Wash bottles Nalgene Type 2403, FEP, with screw cap, ETFE, 125 ml, Screw cap diam. 24 mm	1 pz.
LLG09223232	Wash bottles Nalgene Type 2403, FEP, with screw cap, ETFE, 250 ml, Screw cap diam. 24 mm	1 pz.
LLG09223233	Wash bottles Nalgene Type 2403, FEP, with screw cap, ETFE, 500 ml, Screw cap diam. 28 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

