



Support rods galvaniser steel

Support rods galvaniser steel

Galvanised steel. External diameter 12 mm. With or without M10 thread.



Code	Description	Packaging
LLG04008491	Support rods galvaniser steel, Type With thread, Length 1000 mm	1 pz.
LLG09224110	Support rods galvaniser steel, Type With thread, Length 1000 mm	1 pz
LLG09224130	Support rods galvaniser steel, Type Without thread, Length 1000 mm	1 pz.
LLG09224115	Support rods galvaniser steel, Type With thread, Length 1500 mm	1 pz.
LLG09224135	Support rods galvaniser steel, Type Without thread, Length 1500 mm	1 pz.
LLG04008490	Support rods galvaniser steel, Type With thread, Length 600 mm	1 pz.
LLG09224106	Support rods galvaniser steel, Type With thread, Length 600 mm	1 pz
LLG09224126	Support rods galvaniser steel, Type Without thread, Length 600 mm	1 pz.
LLG06242703	Support rods galvaniser steel, Type With thread, Length 750 mm	1 pz.
LLG09224107	Support rods galvaniser steel, Type With thread, Length 750 mm	1 pz
LLG09224127	Support rods galvaniser steel, Type Without thread, Length 750 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

