

**Culture media bottles, glass, DURAN****DWK Life Sciences****Culture media bottles, glass, DURAN®**

Beaded rim. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable.



Code	Description	Packaging
LLG09010273	Culture media bottles, glass, DURAN®, 5000 ml, Diam. 185 mm, Neck diam. 54 mm, Height 390 mm	1 pz
LLG09010224	Culture media bottles, glass, DURAN®, 100 ml, Diam. 50 mm, Neck diam. 29 mm, Height 115 mm	1 pz
LLG09010254	Culture media bottles, glass, DURAN®, 1000 ml, Diam. 105 mm, Neck diam. 46 mm, Height 237 mm	1 pz
LLG09010244	Culture media bottles, glass, DURAN®, 500 ml, Diam. 83 mm, Neck diam. 42 mm, Height 207 mm	1 pz
LLG09010266	Culture media bottles, glass, DURAN®, 2500 ml, Diam. 150 mm, Neck diam. 50 mm, Height 315 mm	1 pz
LLG09010239	Culture media bottles, glass, DURAN®, 300 ml, Diam. 70 mm, Neck diam. 42 mm, Height 168 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

