

**LLG-Quantitative filter paper, circles****LLG Labware****LLG-Quantitative filter paper, circles**

Made from 100% cellulose.

Ash content 0,007%.

Filtration speed according to DIN 53137.

Capacity of separation according to DIN 53135.

Medium fast, pore size 5 to 8µm.

-for NH_4MgPO_4 , CaC_2O_4 , BaSO_4

-ashless, equivalent red tape



Code	Description	Packaging
LLG06242748	Quantitative filter paper, circles, Diam. 110 mm, Filtration time 50 secs.	100 pz
LLG07970045	Quantitative filter paper, circles, Diam. 125 mm, Filtration time 50 secs.	100 pz
LLG06252765	Quantitative filter paper, circles, Diam. 150 mm, Filtration time 50 secs.	100 pz
LLG06242747	Quantitative filter paper, circles, Diam. 185 mm, Filtration time 50 secs.	100 pz
LLG07970591	Quantitative filter paper, circles, Diam. 240 mm, Filtration time 50 secs.	100 pz.
LLG06256437	Quantitative filter paper, circles, Diam. 55 mm, Filtration time 50 secs.	100 pz
LLG07970571	Quantitative filter paper, circles, Diam. 70 mm, Filtration time 50 secs.	100 pz
LLG07970032	Quantitative filter paper, circles, Diam. 90 mm, Filtration time 50 secs.	100 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

