

**Overhead mixer LLG-uniROTATOR 2****LLG Labware****Overhead mixer LLG-uniROTATOR 2**

Analogue overhead mixer with vertical rotary motion for gentle but effective mixing. Ideal for biological samples. Suitable for cutting and mixing processes of powdered and liquid samples in tubes up to 50 ml. With easily interchangeable clamps that allow a quick cleaning of the mixer.

- Compact design
- Housing easy to clean
- Non-slip, adjustable screw feet

Scope of delivery: Overhead mixer incl. 12 clamps each for 2 ml, 15 ml and 50 ml tubes

Specifications

Type of movement:

Speed: overhead

0 ... 80 min -1

Rotation angle: 360°

Permissible shaking weight (incl. attachment):

Permissible relative moisture:

Permissible ambient temperature: 2 kg

80 %

5 ... 40 °C

Motor rating input: 20 W

Dimensions (W x D x H): 530 x 150 x 190 mm

Power supply: 100 ... 240 V, 50/60 Hz

IP code: IP 21

Warranty: 3 years



Code	Description	Packaging
LLG06263660	uniROTATOR 2 Overhead rotator with EU-plug	1 pz
LLG06263661	uniROTATOR 2 Overhead rotator with UK-plug	1 pz.
LLG06263663	Accessories for Overhead mixer uniROTATOR 2 / uniLOOPMIX 2, Clamp for 15 ml centrifuge tube	12 pz.
LLG06263662	Accessories for Overhead mixer uniROTATOR 2 / uniLOOPMIX 2, Clamp for 2 ml centrifuge tube	12 pz.
LLG06263664	Accessories for Overhead mixer uniROTATOR 2 / uniLOOPMIX 2, Clamp for 50 ml centrifuge tube	12 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

