

**Disc rotator LLG-uniLOOPMIX 2****LLG Labware****Disc rotator LLG-uniLOOPMIX 2**

Analogue Disc rotator for gentle and effective sample mixing. Suitable for all biological samples. With easily interchangeable clamps that allow a **quick cleaning** of the rotator.

For a maximum of 48 x 1.5 ml, 16 x 5 or 15 ml, 8 x 50 ml tubes. Also for parallel use of different sizes, can be equipped individually.

- Compact design
- Housing easy to clean
- Non-slip, adjustable screw feet
- Additional clamps on request

Scope of delivery: Disc rotator incl. 12 clamps each for 2 ml, 15 ml and 50 ml tubes

Specifications

Type of movement:

Speed: rotating

0 ... 80 min -1

Rotation angle: 360°

Permissible shaking weight (incl. attachment): 2 kg

Permissible relative moisture:

Permissible ambient temperature: 80 %

5 ... 40 °C

Motor rating input/output: 20 W

Dimensions (W x D x H): 220 x 300 x 310 mm

Power supply: 100 ... 240 V, 50/60 Hz

IP code: IP 21

Warranty: 3 years



Code	Description	Packaging
LLG06263666	uniLOOPMIX 2 Loop rotator with UK-plug	1 pz.
LLG06263665	uniLOOPMIX 2 Loop rotator with EU-plug	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

