



Support rings with or without bosshead, 18/10 Stainless / zincked

Support rings with or without bosshead, 18/10 Stainless / zincked



Code	Description	Packaging
LLG06238844	Closed ring, with bosshead, Int. diam. 100 mm, Length 70 mm, Material steel zincked	1 pz.
LLG06239804	Closed ring, Int. diam. 100 mm, Length 220 mm, Material steel zincked	1 pz.
LLG09225030	Closed ring, Int. diam. 100 mm, Length 220 mm, Material 18/10 steel	1 pz.
LLG06239822	Closed ring, Int. diam. 130 mm, Length 220 mm, Material steel zincked	1 pz.
LLG07651881	Closed ring, with bosshead, Int. diam. 130 mm, Length 70 mm, Material steel zincked	1 pz.
LLG09225033	Closed ring, Int. diam. 130 mm, Length 220 mm, Material 18/10 steel	1 pz.
LLG06239507	Closed ring, Int. diam. 70 mm, Length 200 mm, Material steel zincked	1 pz.
LLG09225027	Closed ring, Int. diam. 70 mm, Length 200 mm, Material 18/10 steel	1 pz.
LLG09225081	Closed ring, with bosshead, Int. diam. 70 mm, Length 70 mm, Material steel zincked	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

