



Accessories for ETS-D5 / ETS-D6 contact thermometer

IKA

Accessories for ETS-D5 / ETS-D6 contact thermometer



Code	Description	Packaging
LLG09232994	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 44, Bosshead	1 pz.
LLG09233000	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 70, Extension cable, 1 metre	1 pz.
LLG09232996	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 38, Holding rod	1 pz.
LLG09910321	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type PC 1.5 Cable, PC Cable	
LLG09720276	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 11, Replacement power cable	1 pz.
LLG09232998	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 62.51, Replacement sensor, stainless steel	1 pz.
LLG09720257	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 16 V, Stainless steel support rod	1 pz.
LLG09232999	Accessories for ETS-D5 / ETS-D6 contact thermometer, Type H 66.51, Stainless steel sensor, glass sleeved	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

