



Cleaning baskets, stainless steel wire

Cleaning baskets, stainless steel wire

With cover, made of non-rusting, electrochemically polished, stainless steel wire. For washing hardware. Frame, base and cover made of 8 x 8 x 1.2mm spot welded wire mesh.



Code	Description	Packaging
LLG09908130	Cleaning baskets, stainless steel wire, Material Stainless steel, Width 150 mm, Depth 250 mm, Height 150 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

