


**8-step irreversible temperature strips testoterm****Testo****8-step irreversible temperature strips testoterm®**

Strips for 8 temperatures. 5 types cover the measuring range 37 to 260 °C, each strip displays 8 different temperatures. Each strip measures 50 x 18mm.



Code	Description	Packaging
LLG09246004	8-step irreversible temperature strips testoterm®	10 pz.
LLG09246005	8-step irreversible temperature strips testoterm®	10 pz.
LLG09246006	8-step irreversible temperature strips testoterm®	10 pz.
LLG09246002	8-step irreversible temperature strips testoterm®	10 pz.
LLG09246003	8-step irreversible temperature strips testoterm®	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

