



## 5-step irreversible temperature strips Thermax

### 5-step irreversible temperature strips Thermax®

5-step irreversible temperature measuring strips covering the temperature range 29 to 290°C. Each strip measures 35 x 15mm.



| Code        | Description                                             | Packaging |
|-------------|---------------------------------------------------------|-----------|
| LLG09246318 | 5-step irreversible temperature strips Thermax®, Type H | 10 pz     |
| LLG09246314 | 5-step irreversible temperature strips Thermax®, Type D | 10 pz     |
| LLG09246315 | 5-step irreversible temperature strips Thermax®, Type E | 10 pz     |
| LLG09246316 | 5-step irreversible temperature strips Thermax®, Type F | 10 pz     |
| LLG09246317 | 5-step irreversible temperature strips Thermax®, Type G | 10 pz     |
| LLG09246319 | 5-step irreversible temperature strips Thermax®, Type I | 10 pz     |
| LLG09246311 | 5-step irreversible temperature strips Thermax®, Type A | 10 pz     |
| LLG09246312 | 5-step irreversible temperature strips Thermax®, Type B | 10 pz     |
| LLG09246313 | 5-step irreversible temperature strips Thermax®, Type C | 10 pz     |
| LLG09246310 | 5-step irreversible temperature strips Thermax®, Type S | 10 pz     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

