



5-step irreversible temperature strips Thermax

5-step irreversible temperature strips Thermax®

5-step irreversible temperature measuring strips covering the temperature range 29 to 290°C. Each strip measures 35 x 15mm.



Code	Description	Packaging
LLG09246318	5-step irreversible temperature strips Thermax®, Type H	10 pz
LLG09246314	5-step irreversible temperature strips Thermax®, Type D	10 pz
LLG09246315	5-step irreversible temperature strips Thermax®, Type E	10 pz
LLG09246316	5-step irreversible temperature strips Thermax®, Type F	10 pz
LLG09246317	5-step irreversible temperature strips Thermax®, Type G	10 pz
LLG09246319	5-step irreversible temperature strips Thermax®, Type I	10 pz
LLG09246311	5-step irreversible temperature strips Thermax®, Type A	10 pz
LLG09246312	5-step irreversible temperature strips Thermax®, Type B	10 pz
LLG09246313	5-step irreversible temperature strips Thermax®, Type C	10 pz
LLG09246310	5-step irreversible temperature strips Thermax®, Type S	10 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

