

**Electric Bunsen burner BA6101****Electrothermal****Electric Bunsen burner BA6101**

Electrothermal's Electric Bunsen, combines the advantages of a regular gas burner with the clean, easy operation of their Electromantles. The Electric Bunsen is ideal for non-contact heating of test tubes, crucibles, small flasks and beakers, regardless of their shape.

- Conical shaped heating element directs radiant heat to the top cavity
- Top cowl deflects heat away from your hand
- Air circulation from the vented housing keeps the base cool enough to hold during operation
- Corrosion-resistant

**Specifications**

**Power consumption: 400W**

**Dimensions (diam. x H): 120 x 177mm**

**Weight: 0.5kg**

**Temperture range: 800 to 1000°C**



Code	Description	Packaging
LLG04011122	Electric Bunsen burner BA6101, Type BA6101/C, with controller	1 pz.
LLG09018691	Electric Bunsen burner BA6101, Type BA6101, without controller	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

