

**Electric Bunsen burner BA6101****Electrothermal****Electric Bunsen burner BA6101**

Electrothermal's Electric Bunsen, combines the advantages of a regular gas burner with the clean, easy operation of their Electromantles. The Electric Bunsen is ideal for non-contact heating of test tubes, crucibles, small flasks and beakers, regardless of their shape.

- Conical shaped heating element directs radiant heat to the top cavity
- Top cowl deflects heat away from your hand
- Air circulation from the vented housing keeps the base cool enough to hold during operation
- Corrosion-resistant

Specifications

Power consumption: 400W

Dimensions (diam. x H): 120 x 177mm

Weight: 0.5kg

Temperature range: 800 to 1000°C



Code	Description	Packaging
LLG04011122	Electric Bunsen burner BA6101, Type BA6101/C, with controller	1 pz.
LLG09018691	Electric Bunsen burner BA6101, Type BA6101, without controller	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

