



Funnels, DURAN

Funnels, DURAN®

DIN ISO 4798. With short stem. Glass Type I / neutral glass as per USP, EP and JP.
Autoclavable.



Code	Description	Packaging
LLG09251123	Funnels, DURAN®, Diam. 35 mm, For filter paper 45 - 55 diam. mm, Height 60 mm, DIN 4798	1 pz
LLG09251133	Funnels, DURAN®, Diam. 55 mm, For filter paper 70 - 90 diam. mm, Height 95 mm, DIN 4798	1 pz
LLG09251157	Funnels, DURAN®, Diam. 150 mm, For filter paper 240 - 270 diam. mm, Height 265 mm, DIN 4798	1 pz
LLG09251146	Funnels, DURAN®, Diam. 100 mm, For filter paper 150 - 185 diam. mm, Height 180 mm, DIN 4798	1 pz
LLG09251023	Funnels, DURAN®, Diam. 120 mm, For filter paper 185 - 240 diam. mm, Height 210 mm, DIN	1 pz
LLG09251161	Funnels, DURAN®, Diam. 200 mm, For filter paper 320 - 385 diam. mm, Height 325 mm, DIN	1 pz
LLG09251169	Funnels, DURAN®, Diam. 300 mm, For filter paper 500 diam. mm, Height 409 mm, DIN	1 pz
LLG09251128	Funnels, DURAN®, Diam. 45 mm, For filter paper 55 - 70 diam. mm, Height 80 mm, DIN	1 pz
LLG09251138	Funnels, DURAN®, Diam. 70 mm, For filter paper 110 - 125 diam. mm, Height 125 mm, DIN	1 pz
LLG09251141	Funnels, DURAN®, Diam. 80 mm, For filter paper 125 - 150 diam. mm, Height 140 mm, DIN	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

