

**Analytical funnels, DURAN****DWK Life Sciences****Analytical funnels, DURAN®**

For rapid filtration. Glass Type I / neutral glass as per USP, EP and JP.
Autoclavable.



Code	Description	Packaging
LLG09251241	Analytical funnels, DURAN®, Funnel diam. 80 mm, For filter paper 110 to 125 diam. mm, Height 210 mm	1 pz
LLG09251237	Analytical funnels, DURAN®, Funnel diam. 65 mm, For filter paper 70 to 90 diam. mm, Height 200 mm	1 pz
LLG09251248	Analytical funnels, DURAN®, Funnel diam. 110 mm, For filter paper 150 to 185 diam. mm, Height 265 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

