

**Filter funnels, PP, Buchner****Kartell****Filter funnels, PP, Buchner**

PP. Two-piece construction. Top and bottom sections can be separated. Easy to clean. Autoclavable. Suitable for foodstuff.



Code	Description	Packaging
LLG09252334	Filter funnels, PP, Buchner, For filterpaper diam. 90 mm, 390 ml	1 pz
LLG09252332	Filter funnels, PP, Buchner, For filterpaper diam. 70 mm, 180 ml	1 pz
LLG09252338	Filter funnels, PP, Buchner, For filterpaper diam. 160 mm, 2100 ml	1 pz
LLG09252333	Filter funnels, PP, Buchner, For filterpaper diam. 80 mm, 285 ml	1 pz
LLG09252330	Filter funnels, PP, Buchner, For filterpaper diam. 42.5 mm, 40 ml	1 pz
LLG09252339	Filter funnels, PP, Buchner, For filterpaper diam. 240 mm, 6000 ml	1 pz
LLG09252331	Filter funnels, PP, Buchner, For filterpaper diam. 55 mm, 70 ml	1 pz
LLG09252335	Filter funnels, PP, Buchner, For filterpaper diam. 110 mm, 810 ml	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

