

**4-Position-Magnetic Stirrer LLG-uniSTIRRER 1/M4****LLG Labware****4-Position-Magnetic Stirrer LLG-uniSTIRRER 1/M4**

An **Ultra-flat** inductive magnetic stirrer without heating with timer. The modern magnet coil technology consists only of non-moving parts and therefore guarantees a wear-free drive.

For a better mixing of samples it is possible to change the direction of rotation every 30 seconds.

- Max. stirring quantity per stirring position 400 ml
- Timer** (1 to 99 min, infinitely)
- 100 % wear- and maintenance-free inductive drive
- LED-display for an exact adjustment of the rotation speed
- Digital regulated rotation speed
- Alteration of the stirring direction
- High Protection class IP 65
- Stainless steel surface

Scope of supply: magnetic stirrer, 1 mains adapter each US, AUS, UK, 1 stirring rod 25 x 8 mm

Specifications**Stirring positions:**

Distance between Stirring positions:

**Max. Stirring quantity (H₂O): 4
95 mm**

4 x 0.4 l = 1.6 l

Speed range: 15 ... 1200 min⁻¹

Permissible ambient temperature: 5 ... 40 °C

Magnetic bar: max. 25 mm

Display: LED

Heating plate (W x D): 122 x 385 mm


Dimensions (W x D x H): 126 x 457 x 13.5 mm

Weight: 1.9 kg

Power supply: 100 ... 240 V, 50/60 Hz

IP code: IP 65

Warranty: 3 years



Code	Description	Packaging
LLG06263425	uniSTIRRER 1/M4 4-position magnetic stirrer with EU plug	1 pz
LLG06263426	uniSTIRRER 1/M4 4-position magnetic stirrer with UK plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

