

**LLG-Centrifuge Tubes economy, PP****LLG Labware****LLG-Centrifuge Tubes economy, PP**

- With flat PE caps
- With conical bottom (15 and 50 ml) or with self-standing bottom (50 ml)
- Easy-to-read black volume graduations
- With a large frosted writing area
- Engraved graduation at the conical bottom of each tube
- Max. rotate speed up to 12000 x g for conical bottom tubes and 6000 x g for self-standing tubes
- Autoclavable at 121 °C and freezable to -80 °C
- Leak-proof
- Free of heavy metals, pyrogens and endotoxins
- DNase/RNase-free
- Aseptically or sterile



Code	Description	Packaging
LLG04668482	LLG-Centrifuge tubes economy conical, PP, 15 ml, sterile, per bag of 25	500 pz
LLG04668483	LLG-Centrifuge tubes economy conical, PP, 50 ml, sterile, per bag of 25	500 pz
LLG04668486	LLG-Centrifuge tubes economy with skirt, PP, 50 ml, sterile, per bag of 25	500 pz
LLG06270403	LLG-Centrifuge tubes economy conical, PP, 15 ml, aseptic, bulk	500 pz
LLG06270404	LLG-Centrifuge tubes economy conical, PP, 50 ml, aseptic, bulk	500 pz
LLG06270405	LLG-Centrifuge tubes economy with skirt, PP, 50 ml, aseptic, bulk	500 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

