


**Watch glasses, PP****Kartell****Watch glasses, PP**

With base. Resistant to steam sterilisation (121°C). Autoclavable.



Code	Description	Packaging
LLG09263530	Watch glasses, PP, Diam. 100 mm	1 pz
LLG09263532	Watch glasses, PP, Diam. 118.5 mm	1 pz
LLG09263526	Watch glasses, PP, Diam. 60 mm	1 pz
LLG09263528	Watch glasses, PP, Diam. 80 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

