



## Single channel microliter pipettes LLG-proMLP

## LLG Labware

The LLG-proMLP single channel microliter pipettes combine high accuracy and reliability with ergonomic and extremely low force pipetting relieving stress from your hand and thumb. This high precision instrument has been built to maintain accuracy for many years making it the ideal choice for all scientific work requiring exact and reproducible pipetting.

The LLG-proMLP pipettes are equipped with an adjustable eject button allowing right and left handed users to set the button\* according to their preferences. The built in volume locking mechanism and competitive pricing complete these high quality precision pipettes. Available volume ranges from 0.2 µl to 10 ml.

- High quality **Made in France**
- High accuracy, reliability and reproducibility
- Extremely light weight while the robustness still remains
- Perfect ergonomics for highest comfort
- Colour-coded push-button for quick operation
- Extremely low pipetting forces
- **Built-in volume locking** for error-free results
- **Tip ejector button\*, adjustable for right- and left-handed users**
- Completely autoclavable at 121°C
- Including individual calibration certificate (according to EN ISO 8655)

\* not available for 5 ml and 10 ml pipettes



Code	Description	Packaging
LLG06270420	Single channel microliter pipettes LLG-proMLP, Volume 0,2-2 µl	1 pz
LLG06270421	Single channel microliter pipettes LLG-proMLP, Volume 0,5-10 µl	1 pz
LLG06270422	Single channel microliter pipettes LLG-proMLP, Volume 2-20 µl	1 pz
LLG06270423	Single channel microliter pipettes LLG-proMLP, Volume 10-100 µl	1 pz
LLG06270424	Single channel microliter pipettes LLG-proMLP, Volume 20-200 µl	1 pz
LLG06270425	Single channel microliter pipettes LLG-proMLP, Volume 100-1000 µl	1 pz
LLG06270426	Single channel microliter pipettes LLG-proMLP, Volume 500-5000 µl	1 pz
LLG06270427	Single channel microliter pipettes LLG-proMLP, Volume 1000-10000 µl	1 pz

## Accessories

Code	Description	Packaging
LLG06288219	Filters for Single channel microliter pipettes proMLP 5K and 10K	10 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*