

**Viscometers Ubbelohde, glass, kinematic****SI Analytics****Viscometers Ubbelohde, glass, kinematic**

ISO 3105, DIN 51562.

For manual (type 501 ..) or automatic (type 532 ..) measurements.

Constants documented in manufacturer's certificate.



Code	Description	Packaging
LLG07023176	Viscometers Ubbelohde, glass, kinematic, Type 532 00, Capillary i.d. 0.36 mm	
LLG09268100	Viscometers Ubbelohde, glass, kinematic, Type 501 00, Capillary i.d. 0.36 mm	1 pz.
LLG09268101	Viscometers Ubbelohde, glass, kinematic, Type 501 01, Capillary i.d. 0.53 mm	1 pz.
LLG09268201	Viscometers Ubbelohde, glass, kinematic, Type 532 01, Capillary i.d. 0.53 mm	
LLG09268203	Viscometers Ubbelohde, glass, kinematic, Type 532 03, Capillary i.d. 0.47 mm	
LLG09268110	Viscometers Ubbelohde, glass, kinematic, Type 501 10, Capillary i.d. 0.63 mm	1 pz.
LLG09268210	Viscometers Ubbelohde, glass, kinematic, Type 532 10, Capillary i.d. 0.63 mm	1 pz.
LLG09268113	Viscometers Ubbelohde, glass, kinematic, Type 501 13, Capillary i.d. 0.84 mm	1 pz.
LLG09268213	Viscometers Ubbelohde, glass, kinematic, Type 532 13, Capillary i.d. 0.84 mm	1 pz.
LLG09268111	Viscometers Ubbelohde, glass, kinematic, Type 501 11, Capillary i.d. 0.95 mm	1 pz.
LLG09268120	Viscometers Ubbelohde, glass, kinematic, Type 501 20, Capillary i.d. 1.13 mm	1 pz.
LLG09268220	Viscometers Ubbelohde, glass, kinematic, Type 532 20, Capillary i.d. 1.13 mm	1 pz.
LLG07023178	Viscometers Ubbelohde, glass, kinematic, Type 532 21, Capillary i.d. 1.69 mm	
LLG09268121	Viscometers Ubbelohde, glass, kinematic, Type 501 21, Capillary i.d. 1.69 mm	1 pz.
LLG09268130	Viscometers Ubbelohde, glass, kinematic, Type 501 30, Capillary i.d. 2.01 mm	1 pz.
LLG09268230	Viscometers Ubbelohde, glass, kinematic, Type 532 30, Capillary i.d. 2.01 mm	1 pz.
LLG09268133	Viscometers Ubbelohde, glass, kinematic, Type 501 33, Capillary i.d. 2.65 mm	1 pz.
LLG09268233	Viscometers Ubbelohde, glass, kinematic, Type 532 33, Capillary i.d. 2.65 mm	1 pz.
LLG09268131	Viscometers Ubbelohde, glass, kinematic, Type 501 31, Capillary i.d. 3.00 mm	1 pz.
LLG07023179	Viscometers Ubbelohde, glass, kinematic, Type 532 31, Capillary i.d. 3.00 mm	



Code	Description	Packaging
LLG09268140	Viscometers Ubbelohde, glass, kinematic, Type 501 40, Capillary i.d. 3.60 mm	1 pz.
LLG09268240	Viscometers Ubbelohde, glass, kinematic, Type 532 40, Capillary i.d. 3.60 mm	
LLG09268103	Viscometers Ubbelohde, glass, kinematic, Type 501 03, Capillary i.d. 0.47 mm	1 pz.
LLG09268123	Viscometers Ubbelohde, glass, kinematic, Type 501 23, Capillary i.d. 1.50 mm	1 pz.
LLG09268223	Viscometers Ubbelohde, glass, kinematic, Type 532 23, Capillary i.d. 1.50 mm	1 pz.
LLG07023177	Viscometers Ubbelohde, glass, kinematic, Type 532 11, Capillary i.d. 0.95 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

