

**Sampler LiquiSampler, PTFE/FEP or PP****Bürkle****Sampler LiquiSampler, PTFE/FEP or PP**

DIN EN ISO 15528. For liquids. Made of ultra pure and chemically inert PTFE/FEP or PP for contamination-free sampling. Practical push button, single-handed operation. Suitable for point, cross-sectional, and bottom sampling. Ideal for sampling from open and closed drums, kegs, tanks, silos, water sources (ponds, lakes, streams, rivers). Suitable for all commercial kegs and drums with openings of at least 32 mm diameter. All surfaces are non-porous and crevice-free preventing accumulation of dirt. Only round-profile screwthreads, suitable for food hygiene applications, are used. Simple to dismantle and clean.



Code	Description	Packaging
LLG06255307	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 150 ml, Length 600 mm, Diam. 32 mm, Colour highly transparent, Material PTFE/FEP	1 pz.
LLG09303950	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 250 ml, Length 1000 mm, Diam. 32 mm, Colour highly transparent, Material PTFE/FEP	1 pz.
LLG09303951	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 500 ml, Length 2000 mm, Diam. 32 mm, Colour highly transparent, Material PTFE/FEP	1 pz.
LLG06901807	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 150 ml, Length 600 mm, Diam. 32 mm, Colour transparent, Material PP	1 pz.
LLG09303955	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 250 ml, Length 1000 mm, Diam. 32 mm, Colour transparent, Material PP	1 pz.
LLG09303956	Sampler LiquiSampler, PTFE/FEP or PP, Nominal capacity 500 ml, Length 2000 mm, Diam. 32 mm, Colour transparent, Material PP	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

