

**Core Sampler Quick Picker, Stainless steel V4A or PP****Bürkle****Core Sampler Quick Picker, Stainless steel V4A or PP**

Suitable for removing bulk materials from sacks or large bags and depositing them directly into sample bottles. Simple to use and clean. Crevice-free to prevent sample residues causing cross-contamination which could affect subsequent sample results. Complete with cleaning brush and two 250 ml sample bottles.

Lock seal close-it and close-it food & pharma please order separately.



Code	Description	Packaging
LLG06205458	Core Sampler Quick Picker, Stainless steel V4A or PP, 75 ml, Length 50 cm, Imm. depth 300 mm, Material PP clear	1 pz.
LLG09303850	Core Sampler Quick Picker, Stainless steel V4A or PP, 75 ml, Length 50 cm, Imm. depth 300 mm, Material V4A	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

