

**Control seals for sampling close-it****Bürkle****Control seals for sampling close-it**

White/plain or coloured with labelling field . close-it has a sandwich-type structure (aluminium/ paper) and seals the container hermetically (vapour barrier). The rear of the label is coated with a powerful adhesive. This means that close-it sticks firmly to most commercial containers, even to surfaces coated with fine powder, slightly moist or even frozen surfaces where conventional labels or adhesive tape will not stick.

close-it food & pharma ensures aswell the necessary success when resealing sample openings, as it reliably prevents products from escaping after sampling. The adhesive complies with the European guidelines for foodstuffs and FDA guidelines. Just like the original, close-it food & pharma has a barrier layer of aluminium foil, which makes the close-it food & pharma perfectly suited to sampling foodstuffs, pharmaceuticals, cosmetics, food additives, even frozen goods - any applications in which the product must not be contaminated by contact with the adhesive.



Code	Description	Packaging
LLG09303914	Control seals for sampling close-it, Type close-it, Printed, 95 x 95mm, Colour Black	500 pz.
LLG09303984	Control seals for sampling close-it, Type close-it Maxi, Printed, 150 x 150mm, Colour Black	250 pz.
LLG09303999	Control seals for sampling close-it, Type close-it food & pharma, Printed, 95 x 95mm, Colour Black	500 pz.
LLG09304003	Control seals for sampling close-it, Type close-it food & pharma Maxi, Printed, Roll, 150mm x 150mm, Colour black	1 pz.
LLG09303913	Control seals for sampling close-it, Type close-it, Printed, 95 x 95mm, Colour Blue	500 pz.
LLG09303983	Control seals for sampling close-it, Type close-it Maxi, Printed, 150 x 150mm, Colour Blue	250 pz.
LLG09303997	Control seals for sampling close-it, Type close-it food & pharma, Printed, 95 x 95mm, Colour Blue	500 pz.
LLG09304001	Control seals for sampling close-it, Type close-it food & pharma Maxi, Printed, Roll, 150mm x 150mm, Colour blue	1 pz.
LLG09303911	Control seals for sampling close-it, Type close-it, Printed, 95 x 95mm, Colour Green	500 pz.
LLG09303981	Control seals for sampling close-it, Type close-it Maxi, Printed, 150 x 150mm, Colour Green	250 pz.
LLG09303998	Control seals for sampling close-it, Type close-it food & pharma, Printed, 95 x 95mm, Colour Green	500 pz.
LLG09304002	Control seals for sampling close-it, Type close-it food & pharma Maxi, Printed, Roll, 150mm x 150mm, Colour green	250 pz.
LLG06228255	Control seals for sampling close-it, Type close-it food & pharma, Printed, 95 x 95mm, Colour Red	500 pz.
LLG09303910	Control seals for sampling close-it, Type close-it, Printed, 95 x 95mm, Colour Red	500 pz.
LLG09303980	Control seals for sampling close-it, Type close-it Maxi, Printed, 150 x 150mm, Colour Red	250 pz.
LLG09303990	Control seals for sampling close-it, Type close-it food & pharma Maxi, Printed, 150 x 150mm, Colour Red	250 pz.
LLG09303915	Control seals for sampling close-it, Type close-it, Plain, 95 x 95mm, Colour White	500 pz.
LLG09303985	Control seals for sampling close-it, Type close-it Maxi, Plain, 150 x 150mm, Colour White	250 pz.
LLG09303986	Control seals for sampling close-it, Type close-it food & pharma, Plain, 95 x 95mm, Colour White	500 pz.



Code	Description	Packaging
LLG09303991	Control seals for sampling close-it, Type close-it food & pharma Maxi, Plain, 150 x 150mm, Colour White	250 pz.
LLG09303992	Control seals for sampling close-it, Type close-it tape, Roll, 95mm x 50m, Colour White	1 pz.
LLG09303993	Control seals for sampling close-it, Type close-it tape Maxi, Roll, 150mm x 50m, Colour White	1 pz.
LLG09303912	Control seals for sampling close-it, Type close-it, Printed, 95 x 95mm, Colour Yellow	500 pz.
LLG09303982	Control seals for sampling close-it, Type close-it Maxi, Printed, 150 x 150mm, Colour Yellow	250 pz.
LLG09303996	Control seals for sampling close-it, Type close-it food & pharma, Printed, 95 x 95mm, Colour Yellow	500 pz.
LLG09304000	Control seals for sampling close-it, Type close-it food & pharma Maxi, Printed, Roll, 150mm x 150mm, Colour yellow	250 pz.
LLG09303989	Control seals for sampling close-it, Type Roller for close-it Maxi, -, Colour -	1 pz.
LLG09303995	Control seals for sampling close-it, Type close-it tape food & pharma Maxi, Roll, 150mm x 50m, Colour White	1 pz.
LLG09303994	Control seals for sampling close-it, Type close-it tape food & pharma, Roll, 95mm x 50m, Colour White	1 pz.
LLG09303919	Control seals for sampling close-it, Type Roller for close-it, -, Colour -	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

