

**Viscometer holders****SI Analytics****Viscometer holders**

Type 05392. VA steel.

Suitable for all Ubbelohde viscometers without TC sensors. For manual and automatic measurements. Holds the viscometer perpendicularly. Accuracy 1°.

Protects the viscometer from damage.



Code	Description	Packaging
LLG09268790	Viscometer holders, Type 05392	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

