

**Wide mouth jars, PET, with screw cap, PP**

behr

Wide mouth jars, PET, with screw cap, PP

Transparent with white polypropylene cap. Food-safe.
Round, except 1000 ml and 2000 ml which are square.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09103140 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 2000 ml | 10 pz. |
| LLG06242746 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 1000 ml | 10 pz. |
| LLG09103136 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 100 ml | 10 pz. |
| LLG09103148 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 2500 ml | 10 pz. |
| LLG09103149 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 3000 ml | 10 pz. |
| LLG09103137 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 250 ml | 10 pz. |
| LLG09103138 | Wide mouth jars, PET, with screw cap, PP, Nominal capacity 500 ml | 10 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

