

**Disposable spoon spatula LaboPlast / SteriPlast, PS****Bürkle****Disposable spoon spatula LaboPlast® / SteriPlast®, PS**

White. The spoon spatula is a double spatula with a 0.5 ml spoon and a 17 mm spoon blade. The spoon spatula has many different uses in labs, as a sampling spoon, for measuring or transferring small amounts, as a scraper or for mixing and stirring. Specially designed for disposable use, resistant to chemicals and corrosion-proof.

- SteriPlast®: Individually packed, additionally gamma-sterilized
- Production, assembly and packaging according to Cleanroom class 7 (10000)
- According to EU food and FDA guidelines



Code	Description	Packaging
LLG06284543	Disposable spoon spatula LaboPlast® / SteriPlast®, PS, Type LaboPlast®, Length 178 mm, Sterile -	100 pz.
LLG06284544	Disposable spoon spatula LaboPlast® / SteriPlast®, PS, Type SteriPlast®, Length 178 mm, Sterile +	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

