



Control units for mixers MIXdrive

43953

Control units for mixers MIXdrive

MIXcontrol eco: Prize-optimized control unit, speed setting and ON/OFF by turning 2 units.

MIXcontrol 20: Stainless steel housing, digital speed display, 4-step power setting, 20 watts.

MIXcontrol 40: Identical to MIXcontrol 20, but 10-step power setting, 40 watts, 2 sockets for MIXdrive stirres.

Specifications

Display: digital



Code	Description	Packaging
LLG09645970	Control unit MIXcontrol 20 for stirring drive MIXdrive	1 pz
LLG09645971	Control unit MIXcontrol 40 for stirring drive MIXdrive	1 pz
LLG09645969	Control unit MIXcontrol eco for stirring drive MIXdrive	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

