



## Control units for mixers MIXdrive

43953

**Control units for mixers MIXdrive**

**MIXcontrol eco:** Prize-optimized control unit, speed setting and ON/OFF by turning 2 units.

**MIXcontrol 20:** Stainless steel housing, digital speed display, 4-step power setting, 20 watts.

**MIXcontrol 40:** Identical to MIXcontrol 20, but 10-step power setting, 40 watts, 2 sockets for MIXdrive stirres.

**Specifications**

**Display:** digital



Code	Description	Packaging
LLG09645970	Control unit MIXcontrol 20 for stirring drive MIXdrive	1 pz.
LLG09645971	Control unit MIXcontrol 40 for stirring drive MIXdrive	1 pz.
LLG09645969	Control unit MIXcontrol eco for stirring drive MIXdrive	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

