



Accessories for mixer MIXdrive

Accessories for mixer MIXdrive



Code	Description	Packaging
LLG09645972	Accessories for mixer MIXdrive, distriBOX 4 (4-way header)	1 pz.
LLG09645973	Accessories for mixer MIXdrive, distriBOX 8 (8-way header)	1 pz.
LLG09645974	Accessories for mixer MIXdrive, mixBATH (stainless steel bath, electropolished), fits for MIXdrive 6/15/60 (also HT)	1 pz.
LLG06255612	Accessories for mixer MIXdrive, mixBATHcover 15, PP, for mixBath in combination with MIXdrive 15	1 pz.
LLG07639439	Accessories for mixer MIXdrive, mixBATHcover 6, PP, for mixBath in combination with MIXdrive 6	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

