

**Disposable double scoop spatula LaboPlast / SteriPlast, PS****Bürkle****Disposable double scoop spatula LaboPlast® / SteriPlast®, PS**

White. Double scoop spatula with rounded ends. Content 1.0 and 0.25 ml, spatula width 11 mm on one side and 5.7 on the other side. The micro spatula can be used to sample, transfer or process small amounts of chemicals, powders, granulates, pastes, creams or liquids. The sturdy and convenient centre handle facilitates precise work.

- SteriPlast®: Individually packed, additionally gamma-sterilized
- Production, assembly and packaging according to Cleanroom class 7 (10000)
- According to EU food and FDA guidelines



Code	Description	Packaging
LLG06284545	Disposable double scoop spatula LaboPlast® / SteriPlast®, PS, Type LaboPlast®, Length 180 mm, Sterile -	100 pz.
LLG06284546	Disposable double scoop spatula LaboPlast® / SteriPlast®, PS, Type SteriPlast®, Length 180 mm, Sterile +	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

