


**LLG-Storage Boxes, PP****LLG Labware****LLG-Storage Boxes, PP**

* With alphanumeric coding on the edge and the bottom of each cavity. Suitable for storage in the fridge.



Code	Description	Packaging
LLG09405760	Storage Boxes, PP, For ND8, 9, 10, 11, Colour neon-green, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 67 x 67 mm, Array 4 x 4	1 pz.
LLG09405762	Storage Boxes, PP, For ND8, 9, 10, 11, Colour blue, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 67 x 67 mm, Array 4 x 4	1 pz.
LLG09405765	Storage Boxes, PP, For ND8, 9, 10, 11, Colour transparent, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 67 x 67 mm, Array 4 x 4	1 pz.
LLG09405756	Storage Boxes, PP, For ND13, Colour red, Bottle size 4,0 ml, dim. (l x w x h) 130 x 130 x 52* mm, Array 7 x 7	1 pz.
LLG09405750	Storage Boxes, PP, For ND8, 9, 10, 11, Colour blue, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 130 x 130 x 45* mm, Array 9 x 9	1 pz.
LLG09405751	Storage Boxes, PP, For ND8, 9, 10, 11, Colour orange, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 130 x 130 x 45* mm, Array 9 x 9	1 pz.
LLG09405753	Storage Boxes, PP, For ND8, 9, 10, 11, Colour yellow, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 130 x 130 x 45* mm, Array 9 x 9	1 pz.
LLG09405761	Storage Boxes, PP, For ND8, 9, 10, 11, Colour neon-pink, Bottle size 1,5 / 2,0 ml, dim. (l x w x h) 67 x 67 mm, Array 4 x 4	1 pz.
LLG04001528	Storage Boxes, PP, For ND20, Colour blue, Bottle size 5 / 10 / 20 ml, dim. (l x w x h) 130 x 130 x 102 mm, Array 5 x 5	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



