



LLG-PP Screw Seals ND10

LLG Labware

LLG-PP Screw Seals ND10



Code	Description	Packaging
LLG04008234	PP Screw Seals ND10, Caps black, center hole, Septa Silicone white / PTFE beige	100 pz
LLG04008235	PP Screw Seals ND10, Caps black, center hole, Septa Silicone white / PTFE red	100 pz
LLG04008236	PP Screw Seals ND10, Caps black, center hole, Septa PTFE red / Silicone white / PTFE red	100 pz
LLG04008237	PP Screw Seals ND10, Caps black, center hole, Septa Silicone white / PTFE blue, slit	100 pz
LLG06267111	PP Screw Seals ND10, Caps black, center hole, Septa Natural Rubber red-orange / TEF transparent	100 pz.
LLG04008231	PP Screw Seals ND10, Caps black, center hole (no liner), Septa	100 pz
LLG06267112	PP Screw Seals ND10, Caps black, closed, Septa Natural Rubber red-orange / TEF transparent	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

