

**Body Showers/Emergency Showers****B-SAFETY****Body Showers/Emergency Showers**

Special showering heads ensure that large water quantities in gentle streams reaches the injured person without damaged the injured skin.

- The showers meet EN 15154-1:2006
- DIN-DVGW tested and certified
- Marked sign according to DIN is included
- Large dimensioned maintenance-free ball tap for large throughput quantities
- Throughput quantity: approx. 60 litres/min.
- Simple assembly
- Stable structure
- Chemical-resistant PE powder painting
- High speed showering heads made from plastic, with improved spray pattern. Corrosion resistant, calcify- and maintenance free, very durable, even emptying.
- Large, easy-to-use release bar
- Comes with all mounting materials
- Connection 3/4"



Code	Description	Packaging
LLG09733812	Body Showers/Emergency Showers, Inlet Surface mount, Length 565 mm, Attachment Over door wall flange	1 pz.
LLG09733810	Body Showers/Emergency Showers, Inlet Surface mount, Length 625 mm, Attachment Wall flange	1 pz.
LLG09733811	Body Showers/Emergency Showers, Inlet Concealed mount, Length 625 mm, Attachment Wall/60 mm - ¾" pipe	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

