

**Ultrasonic Baths Sonorex Digitec DT ¼ F - Basic equipment****Bandelin****Ultrasonic Baths Sonorex Digitec DT ¼ F - Basic equipment**

Flat ultrasonic baths with fast degassing and for sample preparation.

- Sample preparation in laboratory flasks
- Fast degassing of samples at the push of the button
- Homogenising
- Same power density in all ultrasonic baths
- Reproducibility of results in the different bath sizes
- Uniform sonication of samples irrespective of size and arrangement of the flasks

**Scope of delivery:**

**DT 510 F:** flask holder GL 510 F and 250ml Tickopur R 33

**DT 1028 F:** 2 x flask holder GL 510 F and 250ml Tickopur R 33

Spring clamps EK or flask holders ZF are necessary to fix the laboratory flasks.



Code	Description	Packaging
LLG09877862	Ultrasonic bath DT 510 F SONOREX DIGITEC 4,3l, 560W without heating	1 pz.
LLG09877863	Ultrasonic bath DT 1028 F SONOREX DIGITEC 9.5 l, 1280 W, without heating	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

