



LLG-PP Screw Seals ND24 (UltraBond Seals and ready assembled Seals ND24), PP Screw Caps ND24 (empty) and Septa ND22

LLG-PP Screw Seals ND24 (UltraBond Seals and ready assembled Seals ND24), PP Screw Caps ND24 (empty) and Septa ND22



Code	Description	Packaging
LLG04008295	PP Screw Caps ND24 (empty) and Septa ND22, Caps white, center hole, Septa no liner	100 pz
LLG06233339	Screw-cap PP, white with holesilic/ 1000	1000 pz.
LLG07633064	Screwcap PP, white wiht holeButylr/ 1000	1000 pz.
LLG04008293	PP Screw Caps ND24 (empty) and Septa ND22, Caps Ultrabond*, white, center hole, Septa Silicone white / PTFE beige	100 pz
LLG04008292	PP Screw Caps ND24 (empty) and Septa ND22, Caps Ultrabond*, white, closed top, Septa Silicone white / PTFE beige	100 pz
LLG09003541	PP Screw Caps ND24 (empty) and Septa ND22, Caps white, closed, Septa Silicone white / PTFE beige (mounted)	1000 pz
LLG04008296	PP Screw Caps ND24 (empty) and Septa ND22, Caps white, closed top, Septa no liner	100 pz
LLG04008291	PP Screw Caps ND24 (empty) and Septa ND22, Caps, Septa Silicone white / PTFE beige	100 pz
LLG09003542	PP Screw Caps ND24 (empty) and Septa ND22, Caps, Septa Butyl red / PTFE gray	1000 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

